

BLOOM — Full App Handover for Claude Code

Build the complete Bloom web app from scratch

1. WHAT IS BLOOM

Bloom is the world's first AI-powered IVF companion app — built as a Next.js progressive web app (PWA) that works on iPhone like a native app.

It is NOT a tracker. It is an emotional and clinical companion.

Tagline: "You are not alone in this journey." One-liner: "Because the women who go through IVF are formed under pressure — and they are extraordinary." Competitor: Embie (embie.app) — Bloom differentiates on AI, emotional support, design, and Secret Space.

2. TECH STACK

- **Framework:** Next.js 16 with App Router
 - **Styling:** Tailwind CSS v3 + PostCSS
 - **Auth/Data:** Supabase (auth + database) — or localStorage for demo
 - **AI:** Anthropic API — model `claude-sonnet-4-6` — for Nora AI chat
 - **Deployment:** Vercel
 - **GitHub:** sarahlaeddine-ops/bloom-app
 - **Fonts:** Cormorant Garamond (serif/logo) + DM Sans (body) from Google Fonts
 - **PWA:** Manifest + apple-touch-icon so it installs on iPhone home screen
-

3. BRAND IDENTITY

Logo

- Word: bloom ✦
- Font: Cormorant Garamond, italic, weight 300, letter-spacing 0.12em
- Color: #9B6DC5
- The ✦ is a 4-pointed diamond/star mark

Diamond ✦ meaning

Formed under pressure, still beautiful. 4 pillars: Clarity · Precision · Resilience · Hope

Color Palette

bg:	#FAF7F4	(warm off-white – main background)
card:	#FFFFFF	(white cards)
surface:	#F0EBE8	(secondary surface)
border:	#E8E0DB	(borders and dividers)
accent:	#9B6DC5	(primary brand purple/lilac)
rose:	#E07A8A	(secondary – emotional, feminine)
teal:	#4ABFB0	(tertiary – hopeful, calm)
gold:	#C49A3C	(quaternary – warmth, trust)
text:	#1A1014	(near-black body text)
muted:	#7A6880	(secondary text)
dim:	#C5B8CC	(tertiary/placeholder text)
deep:	#7C3AED	(deep purple for hover states)

App Icon

✦ 4-pointed star in lilac #B388D9 on pink-to-lilac gradient (#FDF0F5 → #EDE0F5), rounded corners 112px radius on 512x512.

4. APP STRUCTURE

```

bloom-app/
├─ app/
│   ├─ globals.css      ← Google Fonts import FIRST, then @tailwind
│   ├─ layout.jsx       ← Google Fonts link tags + PWA meta tags in head
│   └─ page.jsx         ← Entry: splash → auth → onboarding → app shell
├─ lib/
│   └─ store.js         ← Auth + user state (localStorage or Supabase)
├─ components/
│   └─ SplashScreen.jsx

```

```
|   ├── AuthScreen.jsx
|   ├── OnboardingScreen.jsx
|   ├── AppShell.jsx      ← Bottom tab bar with 5 tabs
|   └── screens/
|       ├── HomeScreen.jsx
|       ├── CheckInScreen.jsx
|       ├── NoraScreen.jsx
|       ├── InsightsScreen.jsx
|       └── MoreScreen.jsx ← Contains all sub-screens
└── public/
    ├── apple-touch-icon.png
    └── manifest.json
```

5. NAVIGATION

Bottom tab bar — 5 tabs:

Tab	Icon	Screen
Home	●	HomeScreen
Check-in	✦	CheckInScreen
Nora AI	⌋	NoraScreen
Insights	✱	InsightsScreen
More	...	MoreScreen

Max width: 430px, centered, full height mobile layout.

6. SCREENS — FULL SPECIFICATION

SPLASH SCREEN

- Full screen warm bg #FAF7F4
- Centered: bloom in Cormorant Garamond italic 64px #9B6DC5
- Below: ✦ in 32px #9B6DC5

- Below: thin horizontal divider line
 - Below: "Your IVF companion" in DM Sans 11px uppercase letter-spacing #C5B8CC
 - Below: 3 bouncing dots in accent, rose, teal colors
 - Bottom text: "You are not alone in this journey" in #C5B8CC
 - Shows for 2.5 seconds then fades to Auth
-

AUTH SCREEN

- Logo at top: bloom italic + ✦ + "Your IVF companion"
 - White card with rounded corners
 - Toggle tabs: Sign Up / Sign In
 - Fields: Name (signup only), Email, Password
 - Button: "Create Account" or "Sign In" in #9B6DC5
 - Error state in red
 - Footer: "Your data is private and never sold."
-

ONBOARDING SCREEN (7 steps)

Progress bar at top showing completion %.

Step 1 — Welcome

- Logo + welcome message to user by name
- List of 5 key features with ✦ bullets

Step 2 — Clinic

- "Your clinic" — text input
- Placeholder: "e.g. Emirates Fertility Centre"

Step 3 — Protocol

- "Your protocol" — single select list:
 - Antagonist

- Long Lupron
- Mini IVF
- Natural
- Not sure yet

Step 4 — Phase

- "Your current phase" — single select with icon + description:
 -  Stimulation — Daily injections, monitoring scans
 -  Two Week Wait — After transfer, waiting for beta
 -  Post Retrieval — Eggs retrieved, waiting for embryos
 -  Pre Transfer — Preparing for embryo transfer
 -  Planning — Planning my first or next cycle

Step 5 — Stim Day

- Number selector 1–14 with +/- buttons
- Day selector chips 1–14

Step 6 — Mandatory Free Therapy

- Title: "One free session. Mandatory."
- Explanation: "We made it mandatory because most women who need it would never book it on their own."
- Therapist card: Dr. Sarah Mitchell, Reproductive Psychiatry, 12 years IVF support
- Book button in #9B6DC5

Step 7 — Done

- Summary of all entered data
- "Open Bloom  button

HOME SCREEN

Header:

- Left: bloom + logo italic

- Right: "Day 7 · Sarah" small text

Hero card (purple tint bg):

- Greeting: "Good morning, Sarah ✨"
- Large text: "Stimulation Day 7"
- Sub: "Antagonist Protocol · Emirates Fertility Centre"
- Progress bar showing cycle progress Day 1 → Trigger
- "X days to retrieval est." text

Stats row (3 cards):

- Follicles: 11 (4 mature) — accent color
- E2 pg/mL: 1.8k — rose color
- Next Scan: 8AM Tomorrow — teal color

Follicle Map card:

- Title: "Follicle Map"
- Bubble visualization — each follicle shown as a circle sized by mm
- Right ovary: 18, 16, 15, 14, 12, 11mm — accent color border
- Left ovary: 17, 16, 14, 13, 12mm — teal color border
- Mature (≥ 16 mm) filled, others empty
- Legend: Right ovary (purple) / Left ovary (teal)

Medications Today card:

- Gonal-F 225 IU · 9:00 PM · Done ✓ (purple)
- Cetrotide 0.25mg · 8:00 AM · Done ✓ (rose)
- Progynova 2mg · 8AM/8PM · Pending (gold)
- Folic Acid 400mcg · 8:00 AM · Done ✓ (teal)

CHECK-IN SCREEN

Header: "How are you today?" + "Takes about 1 minute"

Mood selector (5 options with emoji):

- 🥰 Hard (rose)
- 😊 Low (gold)
- 😐 Okay (muted)
- 🌱 Hopeful (accent)
- 🌸 Good (teal)

Anxiety level slider (1–5):

- Filled circles 1–5, rose color

Hopefulness slider (1–5):

- Filled circles 1–5, accent color

Symptoms chips (multi-select): Bloating, Cramping, Headache, Nausea, Breast tenderness, Hot flashes, Fatigue, Injection site pain, Back pain, Mood swings, Spotting, Insomnia

OHSS Weight tracker (amber card):

- Warning: "With 11 follicles, track daily. Alert if +2kg in 24hrs."
- Number input for weight in kg

Journal note:

- Textarea: "How are you really feeling today?"

Save button: "Save today check-in" in accent

Done state: + "Check-in saved" success screen

NORA AI SCREEN

Header:

- Nora avatar (accent circle with +)
- "Nora" bold
- "Online · Knows your cycle" in teal

Chat interface:

- User messages: right-aligned, accent background, white text
- Nora messages: left-aligned, white card, dark text
- Nora avatar on left of each message
- Typing indicator: "..." animated

Suggested questions (shown when no messages):

- "What does my E2 mean?"
- "Am I at risk for OHSS?"
- "When is my trigger shot?"
- "I am scared about retrieval"

System prompt for Nora:

You are Nora, the AI companion inside Bloom – a dedicated IVF tracking app.
 The patient is on Stimulation Day 7, Antagonist Protocol.
 Follicles: Right 18,16,15,14,12,11mm Left 17,16,14,13,12mm. E2=1840.
 Be warm, concise, clinically accurate.
 Always remind the user to confirm with their clinic.
 Keep responses under 4 sentences.

API call:

```
fetch("https://api.anthropic.com/v1/messages", {
  method: "POST",
  headers: { "Content-Type": "application/json" },
  body: JSON.stringify({
    model: "claude-sonnet-4-6",
    max_tokens: 500,
    system: SYSTEM_PROMPT,
    messages: conversationHistory
  })
})
```

Input bar: Fixed at bottom, textarea + send button (→)

INSIGHTS SCREEN

Header: "Insights" + "Personalized to Stimulation Day 7"

Category filter (horizontal scroll):

ID	Label	Icon	Color
all	All	✦	#9B6DC5
mental	Mental Health	🕒	#8B7AC5
nutrition	Nutrition	🌱	#C49A3C
science	IVF Science	💎	#9B6DC5
intimacy	Intimacy	💍	#E07A8A
movement	Movement	🌀	#4ABFB0

Active category: solid color background, white text/icon. Inactive: white card, colored icon.

Most popular (horizontal scroll, large cards 200×240px): Cards with colored bg, video badge if video, category tag, title.

For you today (2-column grid): Same card style, smaller.

Article detail view:

- Back button
- Colored hero area with category tag and title
- Body text
- "Got it" button in category color

MORE SCREEN

Profile card at top:

- Avatar circle (initial letter)
- Name + email
- "Edit →" link

Feature grid (2 columns): Each card has:

- Colored icon box with mark

- Bold label
- Small description

ID	Mark	Label	Color	Description
report	↓	My Cycle Report	#9B6DC5	Download and share with clinic
medications	⊙	Medications	#9B6DC5	Injection tracker and log
appointments	🕒	Appointments	#E07A8A	Scans, retrieval, transfer
charts	↗	Charts & Trends	#E07A8A	Hormone trends, follicle progress
tww	🕒	Two Week Wait	#C49A3C	Countdown and daily science
therapy	◇	Therapy & Coaching	#4ABFB0	Book IVF-specialist therapists
videos	▶	Wellbeing Videos	#4ABFB0	Movement, breathwork, meditation
community	⊙	Community	#9B6DC5	Anonymous rooms by IVF phase
failed	◆	After a Failed Cycle	#5BADD4	Grief support and next steps
partner	🌙	Partner Space	#FDBA74	Invite and connect your partner
pregnant	✦	Pregnancy Journey	#E07A8A	Week-by-week pregnancy guide
secret	🔒	Secret Space	#8B7AC5	Private journal, only you can see
upgrade	✦	Upgrade to Bloom+	#9B6DC5	Unlock all features

At a glance card at bottom: Shows Name, Phase, Protocol, Follicles, E2 today

MORE — SUB SCREENS

MEDICATIONS SCREEN

3 tabs: Today / History / Schedule

Today tab:

- Done all banner (teal) if all taken
- Injections section: Gonal-F, Cetrotide
- Oral medications section: Progynova, Folic Acid
- Each med card shows: name, dose, time, next dose, Done/Log badge
- Injection tips card (gold border)

Log dose modal (bottom sheet):

- Taken / Missed toggle
 - Injection site selector (Left belly, Right belly, Left thigh, Right thigh) — only for injections
 - Note textarea
 - Cancel / Save buttons
-

APPOINTMENTS SCREEN

Hero card:

- "Next appointment" label
- "Monitoring Scan" title
- "Tomorrow at 8:00 AM · Emirates Fertility Centre"
- 3-box countdown: days / hours / min
- "How to prepare" section (purple bg, 4 tips)

All upcoming list:

- Monitoring Scan — Tomorrow — 1d (accent)
- Trigger Shot Timing — In 2 days — 2d (rose)
- Egg Retrieval — In ~4 days — 4d (gold)
- Embryo Transfer — In ~9 days — 9d (teal)

- Beta HCG Test — In ~23 days — 23d (accent)

Tap to expand: shows doctor name.

WELLBEING VIDEOS SCREEN

Mood filter (horizontal scroll): I feel anxious, I feel scared, I feel hopeful, I need rest, I can't cope, I need peace, I feel bloated, I am grieving, I want to help, I feel nervous

Category filter: Gentle Movement (teal), Breathwork (accent), Meditation (lavender), Nutrition (gold), Emotional Support (rose)

Video cards (list):

- Colored thumbnail box with category mark
- Phase badge (Stimulation, TWW, Any phase, etc.)
- Title
- Duration
- Heart save button

10 videos covering yoga, breathwork, meditation, nutrition, emotional support.

Video player view:

- Large colored thumbnail
- Play button overlay
- Title, phase, duration
- "Opens on YouTube" note

COMMUNITY SCREEN

Rooms list:

Room	Color	Members
Stimulation	#9B6DC5	124
Two Week Wait	#E07A8A	89

After Retrieval	#C49A3C	56
Failed Cycle	#5BADD4	43
IVF Success	#4ABFB0	201
General	#7A6880	312

Anonymous banner: "Anonymous by default — only your flower name is shown"

Inside a room:

- Messages with flower names (Anonymous Sunflower, Anonymous Rose, etc.)
 - Like button (♡ / ♥)
 - Prompt chips: "How are you feeling today?", "What helped you most this week?", etc.
 - Input bar at bottom with send button
-

CYCLE REPORT SCREEN

Full clinical report with:

- Patient info: name, clinic, doctor, protocol, cycle start, stim day
 - Hormone table: Day / Date / E2 / LH / P4 (across all scan days)
 - Follicle map: right and left ovary, each follicle sized, mature flagged
 - Medications list: name, dose, frequency, type, start day
 - Daily check-in history: date, mood, anxiety, hope, symptoms
 - Disclaimer text
 - "Share with my clinic →" button (uses Web Share API)
 - "Copy as text" button
-

SECRET SPACE SCREEN

Locked state:

-  icon

- "Secret Space"
- "A private space for the feelings you cannot say out loud."
- PIN input (password type)
- Enter button (any 4+ digit PIN unlocks for demo)

Unlocked state:

- "Secret Space — Only visible to you" banner
 - Prompt chips: "I resent friends who got pregnant easily", "I Googled success rates at 3am again", "I am terrified this will never work"
 - Chat interface (same as Nora but without API — simple empathetic responses)
 - Fixed input bar at bottom
-

PARTNER SPACE SCREEN

Unconnected state:

- "Bloom for Partners" header
- "You are always in control" privacy message
- 4 feature rows: Her cycle phase, Upcoming appointments, How to support her today, IVF explained
- Pairing code display (e.g. BL-A4X2) in large letters
- "Share via WhatsApp" button (green)
- "Copy code" button

Connected state:

- "How to support her today" — 3 daily tip cards
 - "What your partner sees" — 4 item list
 - FAQ accordion (4 questions)
 - "Disconnect partner" button (rose)
-

UPGRADE SCREEN (Paywall)

Monthly/Annual toggle:

- Annual shows "Save 46%" badge in teal

Plan cards:

Bloom+ — \$19.99/month / \$129/year Color: #9B6DC5 Features: Unlimited Nora AI, Full Insights library, Wellbeing Videos, Community chat, Two Week Wait mode, Cycle Report download, Partner Space, Check-in history and trends

Bloom Pro — \$29.99/month / \$199/year Color: #E07A8A Badge: "Most popular" Features: Everything in Bloom+, Secret Space, Unlimited therapy booking, Priority Nora responses, Multi-cycle tracking, Custom medication schedules, Early access to new features, Bloom Pro badge in Community

CTA button: "Start Bloom Pro — \$199/yr" in selected plan color **Trial note:** "7-day free trial · Cancel anytime" **Free tier list:** basic tracking, 5 Nora messages/day, 3 articles/week, onboarding + free therapy **"Continue with free plan"** link at bottom

PROFILE SCREEN

- Avatar circle (initial letter, accent bg)
 - Name input
 - Clinic input
 - Protocol selector chips (Antagonist, Long Lupron, Mini IVF)
 - "Save Changes" button → turns teal "Saved ✓" on success
 - "Sign Out" button (red border)
-

7. DATA MODEL (localStorage / Supabase)

```
user = {
  id: string,
  name: string,
  email: string,
  password: string,          // hashed in Supabase
  clinic: string,            // "Emirates Fertility Centre"
  protocol: string,          // "Antagonist"
  phase: string,             // "stimulation"
```

```
stimDay: number,          // 7
follicles: number,        // 11
e2: number,               // 1840
onboarded: boolean,
}
```

8. LANDING PAGE (separate from app)

One-page, no scroll on desktop, centered content.

Content:

1. bloom ✦ — Cormorant Garamond italic 80px #9B6DC5
2. "Your IVF companion" — 11px uppercase #C5B8CC
3. Headline — 52px Cormorant Garamond: "The hardest journey / you will ever take. / *Not alone.*" (italic purple)
4. Sub — "An AI companion built for women going through IVF. Launching soon."
5. Email form — input + "Join the waitlist" button
6. "2,847 women already waiting" — number in Cormorant 22px accent
7. Footer — "© 2025 Bloom · hello@bloomivf.app"

Success state: "You are on the list. ✦"

Background: #FDFAF7 with subtle radial gradient orbs (purple top-left, rose bottom-right)

9. PWA SETUP

public/manifest.json:

```
{
  "name": "Bloom",
  "short_name": "Bloom",
  "description": "Your IVF companion",
  "start_url": "/",
  "display": "standalone",
  "background_color": "#FAF7F4",
  "theme_color": "#9B6DC5",
  "icons": [
```

```
    { "src": "/apple-touch-icon.png", "sizes": "512x512", "type": "image/png" }  
  ]  
}
```

layout.jsx head tags:

```
<link rel="apple-touch-icon" href="/apple-touch-icon.png" />  
<link rel="manifest" href="/manifest.json" />  
<meta name="apple-mobile-web-app-capable" content="yes" />  
<meta name="apple-mobile-web-app-status-bar-style" content="default" />  
<meta name="apple-mobile-web-app-title" content="Bloom" />  
<meta name="theme-color" content="#FAF7F4" />
```

10. TAILWIND CONFIG

```
module.exports = {  
  content: ["../app/**/*.{js,jsx}", "../components/**/*.{js,jsx}"],  
  theme: {  
    extend: {  
      colors: {  
        bloom: {  
          bg: "#FAF7F4",  
          card: "#FFFFFF",  
          surface: "#F0EBE8",  
          border: "#E8E0DB",  
          accent: "#9B6DC5",  
          deep: "#7C3AED",  
          rose: "#E07A8A",  
          teal: "#4ABFB0",  
          gold: "#C49A3C",  
          text: "#1A1014",  
          muted: "#7A6880",  
          dim: "#C5B8CC",  
        }  
      }  
    }  
  },  
  plugins: []  
}
```

11. GLOBALS.CSS

```
@import url("https://fonts.googleapis.com/css2?family=Cormorant+Garamond:ital,wght@1,600,700,800,900;family=DM+Sans:wght@400,500,700,900");

@tailwind base;
@tailwind components;
@tailwind utilities;

* { box-sizing: border-box; margin: 0; padding: 0; }

body {
  font-family: "DM Sans", sans-serif;
  background: #FAF7F4;
  color: #1A1014;
  max-width: 430px;
  margin: 0 auto;
  min-height: 100vh;
}

::-webkit-scrollbar { display: none; }
```

12. IMPORTANT RULES FOR BUILDING

1. Use `"use client"` at top of every component
2. Use `var` not `const` or `let` in logic outside JSX
3. No `gap:` in inline styles — use `marginRight` or `marginBottom`
4. All screens max-width 430px centered — mobile-first
5. Bottom tab bar fixed at bottom, z-index 50
6. Sub-screens open by setting `active` state in `MoreScreen` — no routing needed
7. Back button: `← Back` in accent color, top-left of each sub-screen
8. Google Fonts `@import` MUST be first line of `globals.css` before `@tailwind`
9. Tailwind v3 only (not v4)
10. No external npm packages except `@supabase/supabase-js` and `lucide-react`

13. GITHUB REPOS

- App: `sarahalaeddine-ops/bloom-app`
- Landing page: `sarahalaeddine-ops/bloom-ivf`

Deploy both to Vercel. Auto-deploys on git push to master.

14. SEED RAISE CONTEXT

- Target: \$3M seed round
- Pricing: Free / Bloom+ \$19.99/mo / Bloom Pro \$29.99/mo
- Target conversion: 5% free-to-paid
- 10,000 paying users = \$220,000 MRR
- Waitlist goal: 500 signups before pitch
- Competitor: Embie — 48K users, clinical tracker, no AI, no emotional layer